VWR® Microbiological Air Sampler SAS Super ISO

INSTRUCTION MANUAL

European catalogue numbers:

710-0871

710-0872

710-0869

710-0870

Version: 1.2
Issued: July 2016

Legal address of manufacturer

Europe

VWR International Srl

Via S. Giusto 85

20153 Milano

Italy

Tel.: + 39 0233 203 11

<http://it.vwr.com>

Country of origin Italy

Table of contents

Warning	5
Safety information	5
References	5
Technical features	6
Intended use	6
Principle	6
The basic idea	6
Models available	7
Brief instructions	7
Operating instructions for SAS Super ISO 100 and 180	7
Visual display presentation	8
Installation	8
The practical use of contact plates	8
List of menus and utility sub menus	9
Preliminary inspection	10
Holder adjustment	10
Petri head adapter (optional)	10
Filling 90 mm Petri dishes	10
Tripod installation (optional)	10
Function	11
Software functions	11
To start with the same air volume as the previous sample	11
“SELECT VOLUME” function (to start with one of the eight memorised volume)	12
“PROGRAM VOLUME” function (modification of the value of a stored volume)	13
Using a “SET DELAY” programme to delay instrument start	14
“UTILITY MODE” function	15
“SET TIME”	15
“SITE”	16
“IDENTIFY”	17

“LANGUAGE”	18
“DISPLAY RECORD”	18
“CLEAR RECORD”	19
“REMOTE CONTROL ”	19
“CALIBRATION”	19
“AUTOSWITCH”	20
“MULTI MODE” function	20
“CONFIGURATION” function.....	22
“POWER OFF” function	26
Operations	26
Adjustable plate holders	26
Troubleshooting	27
Accessories	27
Technical service	29
Warranty	29
Compliance with local laws and regulations	29
Equipment disposal	30

Warning

This document is the property of VWR International S.r.l. - Milan - Italy

It may not be duplicated or distributed without the owner's authorisation.

Patent Pending

Manufactured by VWR International S.r.l. - Milan, Italy

Safety information

Please use this device only for the purposes indicated.

The device must be correctly used according to this instructions manual before starting any operation.

Any electrical cables need to be replaced immediately when damaged; never use a damaged or worn electrical cable.

Always disconnect the charger before:

- Repairing or maintenance; these operations must be carried out by qualified staff
- Cleaning the unit

Use original spare parts and accessories for any replacement.

Do not use this device in the presence of explosive gas.

Please follow the guidelines below and read this manual in its entirety to ensure safe operation of the unit.

Be aware that the voltage and frequency of the electrical system are compatible with the requirements of the battery charger.

Never use a non OEM charger to charge the air sampler.
Use of an improper charger may damage the unit.

Carefully read this manual before operating your instrument.

References

FDA - 1987 Guideline on Sterile Drug Products produced by Aseptic Process

ACGIH - Guideline for the Assessment of Bioaerosol in the Indoor Environment

ASTM - Draft Protocol - Committee D22.05.06

USP 23-NF 18 8th Supplement 1116 (May 1998) - Microbiological Evaluation of Clean Rooms and other Controlled Environments

EU Guide for GMP - Manufacture of Sterile Medicinal Products Control of Medicines and Inspection

CEN/TC 243 Norms for Clean Room Technology

Technical features



Feature	SAS ISO 100 Contact	SAS ISO 180 Contact	SAS ISO 100 Petri	SAS ISO 180 Petri
Cat. No.	710-0871	710-0872	710-0869	710-0870
Air flow	100 l/min	180 l/min	100 l/min	180 l/min
Use with	55 mm contact plates	55 mm contact plates	90 mm Petri dishes	90 mm Petri dishes
Portable	Yes	Yes	Yes	Yes
Battery	Rechargeable	Rechargeable	Rechargeable	Rechargeable
Battery life	70 000 litres	40 000 litres	70 000 litres	40 000 litres
Speed sensor	Yes	Yes	Yes	Yes

Intended use

Principle

The Surface Air System (SAS) encompasses several models which use the same principle.

Air is aspirated at a fixed speed for a variable time through a cover which has been machined with a series of specially designed small holes. The resulting laminar airflow is directed onto the agar surface of a “contact plate” (or Petri dish) containing media suitable for the microbiological examination to be performed. When the pre-set sampling cycle is completed, the plate is removed and incubated. The organisms are then visible to the naked eye and can be counted in order to assess the level of contamination.

The basic idea

The major points of the Surface Air System (SAS) are:

1. To use a simple and inexpensive “contact plate” (“SURFAIR PLATE”, “RODAC”) for surface, hands or air control. These plates are very well known and easily available and can be purchased ready poured with different media.
2. To sample a known volume of air for a variable time to provide a large range of sampling volumes.
3. To aspirate air in a laminar flow pattern with sufficient velocity to impact organisms on an agar surface.
4. To accumulate data on the level of hygiene in each environment so that fluctuations can be monitored.

5. To take advantage of advanced electronics for more reliable results in different operating conditions.
6. To have the flexibility to choose between 55 mm contact plates or 90 mm standard disposable Petri dishes.
7. To apply cGLP and cGMP to air sampling operations.
8. To organise sequential sampling to obtain a more representative sample under actual operating conditions.

Models available

SAS SUPER ISO 100 and SAS SUPER ISO 180: Two instruments for two different applications.

The two air samplers have the same performances with the only difference being in the rate of air aspirated:

SAS SUPER ISO 100 = 100 litres of air per minute.

SAS SUPER ISO 180 = 180 litres of air per minute.

The SAS SUPER ISO 180 is appropriate for use in cleanrooms and other applications demanding fast sampling times. In cleanrooms, for example, it is important to test larger volumes of air, because microbial air contamination is very low. The SAS SUPER ISO 180 reduces the time required to obtain a sample.


Both SAS SUPER ISO 100 and SAS SUPER ISO 180 are identified as "SAS SUPER ISO".

Brief instructions

Brief operating instruction for SAS Super ISO 100 and 180

The airflow can be seen on a digital display at the fourth automatic succession presentation each time the instrument is switched on.

1. Touch the screen for at least two seconds to switch the instrument on.

2. Press START to sample the same air volume as the last sampling cycle.
3. To change the volume of air, press MENU, to reach "select volume", then press ENTER and use the arrows to choose the new desired volume of air.
4. Refer to instruction manual to change settings.
5. Press ENTER to confirm selection.
6. Press MENU and use the arrows for menu selection to modify other pre-set parameters.
7. Press ENTER to reach the sub menu of the chosen parameter.
8. Refer to the instruction manual to change settings.
9. Press STOP each time you need to end an action. The unit will then go back to its initial configuration.
10. Press MENU, and then use the arrows to reach POWER OFF, and then ENTER to switch the unit off.


Visual display

(Touch the screen for at least two seconds).

Each time the instrument is turned on, a five second visual display presentation will appear displaying the following information.


N.B: In the following display, "XX" indicates numerical figures.

This information appears in automatic succession.

Installation

The practical use of contact plates

1. Remove the aspirating head.
2. Insert an identified, closed and prepared plate and remove the plate lid.

3. Replace the aspirating head.
4. Select required volume and start the unit. The airflow is directed into the agar surface of the plate.
5. At the end of the cycle, remove the aspirating head.
6. Close and remove the plate.
7. Incubate.
8. Count the colonies, record the results on the microbiological air sampling report and read the results.

List of menus and utility sub menus

START	Press to repeat the same volume of prior sampling
MENU	Press to enter the following sub menus
SELECT VOLUME	The volume of aspirated air can be modified according to 8 programmed values
PROGRAM VOLUME	The 8 programmed values can be modified
SET DELAY	The air sampler can be programmed to start sampling after a programmable delay (1, 2, 3, 4, 5, 6, 7, 8, 9, 10, 20, 30, 40, 50, 60 minutes)
MULTIMODE	The air sampler can be programmed to extend the total sampling time using "sequential interval time" sampling by the following sub menus:
	INTERVAL TIME Can be programmed from 1 to 59 min
	NUMBERS OF RUNS Can be programmed from 1 to 19 air aspirations
	SINGLE RUN VOLUME Can be programmed from 1 to 1999 l
	TOTAL RUNS Shows the total volume of aspirated air
CONFIGURATION	You can use the SAS with eight programmable configurations
UTILITY MODE	
	SET TIME To programme date and time
	SITE To identify the sampling point
	IDENTIFY To identify the operator's name
	LANGUAGE To display text in different languages
	DISPLAY RECORD To show the recorded sampling data
	CLEAR RECORD To clear the recorded sampling data
	REMOTE CONTROL To connect to a PC with an optional device
	CALIBRATION

	To calibrate the volume of aspirated air
	AUTOSWITCH-OFF To disable the autoswitch off
POWER OFF	The air sampler can be switched off by pressing the ENTER button


Preliminary inspection

The apparatus is subject to specific working tests before shipping and it is carefully packed to avoid possible damages during transportation. However, a visible check should be carried out as soon as possible to determine any transit damage. This must be reported immediately. The following procedures should be followed to check that the unit is working properly.

The battery pack of the instrument must be charged for at least 5 hours before testing SAS SUPER ISO 100 and SAS SUPER ISO 180.

After the screen is touched for at least two seconds, an automatic visual presentation appears (see visual display presentation paragraph).

Press START on the screen to start air sampling with the last cycle. The motor will run until the figures in the brackets reach the same figure as displayed to the right and then it will stop.

Holder adjustment

Contact plate and the Petri dish holders can be adjusted (using a hex key) in case the available plates are slightly different in diameter from the standard 55 mm contact plate or 90 mm Petri dish.

Petri head adapter (optional)

An adapter is available in aluminium or stainless steel which allows the use of standard 90 mm Petri dishes with SAS for contact plates.

Filling 90 mm Petri dishes

We recommend filling standard 90 mm disposable Petri dishes with no more than 18 - 20 ml of agar to avoid the medium touching the inside surface of the aspirating head.

Tripod installation (optional)

The SAS SUPER ISO 100 and SAS SUPER ISO 180 can be fixed to a table/floor tripod or a wall support. The screw thread connection is located underneath the unit, between the two front feet.

Functions

Software function

Touch the screen for at least two seconds to switch the unit on.

The air sampler is supplied with eight programmable air volumes ("Select volume") stored in the memory and eight sampler configuration ("Configuration") selectable by the operator.


The selection of the volume is made by pressing the "UP" or "DOWN" arrows when the programme in the relevant sub menu.

The following volumes of air are suggested:


- Contaminated areas (communities, processing rooms, etc) 10 - 200 litres of air
- Normal areas (laboratory benches, houses, etc) 200 - 500 litres of air
- Sterile or high risk areas (cleanrooms, operating theatres, etc) 500 - 1000 litres of air

To start with the same air volume as the previous sample

Switch on the unit by touching the screen for at least two seconds. After the visual display presentation you will see the following display with the last aspirated volume of air:


Press START on the screen to start air sampling with the last cycle. The motor will run until the figures in the brackets reach the same figure as displayed to the right and then it will stop.


“SELECT VOLUME” function (to start with one of the eight memorised volumes)


Switch the unit on by touching the screen for at least two seconds. After the visual display presentation you will see the following display with the last aspirated volume of air.

Press MENU and select the SELECT VOLUME function by pressing the “UP” or “DOWN” arrows.

Press ENTER and select a volume by pressing the “UP” or “DOWN” arrows (for example 1250 litres)

Press ENTER to select the volume and the display will show the selected volume.

08/31/06 11:31

1250 Lt. DELAY 00 MIN

START **MENU**

Press START to aspirate the selected volume of air.

“PROGRAM VOLUME” function (modification of the value of a stored volume)

Using this procedure, it is possible to memorise up to eight different volumes (from 1 to 1999 litres of air).

Switch the unit on by touching the screen for at least two seconds.

Press MENU and select the PROGRAM VOLUME function by pressing the “UP” or “DOWN” arrows.

MENU'

PROGRAM VOLUME

CLEAR **ENTER**

Press ENTER to modify the volumes.

SELECT VOLUME

[1250] Litres

CLEAR **ENTER**

Select a volume by “^” pressing the “UP” or “DOWN” arrows and press ENTER.

MODIFY RUN

[1250]
^

CLEAR **ENTER**

Change the digit indicated by pressing the “UP” or “DOWN” arrows and presses ENTER to go to the next figure.

Repeat until you reach the last digit and then confirm by pressing ENTER.


Note: The maximum air volume is 1999 litres.

Using a “SET DELAY” programme to delay instrument start


This function is used to set a delay time before the sampling cycle.

Switch the unit on by touching the screen for at least two seconds.

Press MENU and select SET DELAY function by pressing the “UP” or “DOWN” arrows.

Press ENTER to modify the delay time.


Select the desired delay time by pressing the “UP” or “DOWN” arrows.


The available delay times are 0 – 1 – 2 – 3 – 4 – 5 – 6 – 7 – 8 – 9 – 10 – 20 – 30 – 40 – 50 – 60 minutes.

Press ENTER to confirm.

When the delay is set (for example 5 minutes), you will see the delay time on the main display.

The message “DELAY” will then flash until the selected time is reached. The display flashes to confirm delayed start has been selected.


“UTILITY MODE” function

In this programme the following SUB-MENUS are listed: Set time, Site, Identify, Language, Display Record, Clear Record, Remote Control, Calibration and Autoswitch.

Always start from UTILITY MODE to reach one of these SUB-MENUS.

Switch on the unit by touching the screen for at least two seconds.


Press MENU and select UTILITY MODE function by pressing the “UP” or “DOWN” arrows.

Press ENTER to enter the sub menu.

“SET TIME”

This option is used to programme day, month, year and time of the day.

Press ENTER to modify the time.

MM/DD/YY HH:mm			
08/31/06		11:33	
^^			
CLEAR	↓	↑	ENTER

Change the number indicated by “^^” by pressing the “UP” or “DOWN” arrows and press ENTER to go to the next number.

MM/DD/YY HH:mm			
08/31/06		11:33	
^^			
CLEAR	↓	↑	ENTER

Repeat until you reach the last number, confirm by pressing ENTER.

“SITE”

This option is used to memorise the name of the site in the file “DISPLAY RECORD”. Site identification should be changed for samples taken at different sites especially if the results are to be printed.

Select SITE function from UTILITY MODE programme.

UTILITY MODE			
SITE			
CLEAR	↓	↑	ENTER

Press ENTER to modify the site.

SITE			
[ABCDEFGH]			
^			
CLEAR	↓	↑	ENTER

Change the figure indicated by “^” by pressing the “UP” or “DOWN” arrows and press ENTER to go to the next figure.

SITE			
[LABOEF GH]			
^			
CLEAR	↓	↑	ENTER

Repeat until you reach the last figure and then confirm by pressing ENTER.

“IDENTIFY”

This option is used to identify the operator. This should be changed if different operators use the sampler and especially if the date is to be printed. The data is recorded in the file “DISPLAY RECORD”.

Select IDENTIFY function from UTILITY MODE programme.

UTILITY MODE			
IDENTIFY			
CLEAR	↓	↑	ENTER

Press ENTER to modify the identity.

IDENTIFY			
[12345678]			
^			
CLEAR	↓	↑	ENTER

Change the figure indicated by “^” by pressing the “UP” or “DOWN” arrows and press ENTER to go to the next figure.


IDENTIFY			
[34345678]			
^			
CLEAR	↓	↑	ENTER

Repeat until you reach the last figure and then confirm by pressing ENTER.


“LANGUAGE”

Menu text can be selected choosing from different languages.

Select LANGUAGE function from UTILITY MODE programme.

Press ENTER to modify the language.


Change the language by pressing the “UP” or “DOWN” arrows and press ENTER to select the desired language; the display will then show the selected language.

Press CLEAR to exit this function.


“DISPLAY RECORD”

The last 300 samples are memorised in the file “DISPLAY RECORD”. Each sample is identified in chronological date order and shows the date, time, operator, site and volume of air sampled (only after a completed run cycle).

Select DISPLAY RECORD function from UTILITY MODE programme.

Press ENTER to display the following data.


You can view all completed cycles by pressing the “UP” or “DOWN” arrows. The symbol “*” is displayed near the volume aspirated if the volume was aspirated with the multimode function.

Press CLEAR to exit this function.


“CLEAR RECORD”

This option is used to delete all the data memorised in the “DISPLAY RECORD”. Before starting this procedure, be sure that all stored data are not required or have been downloaded.

Select CLEAR RECORD function from UTILITY MODE programme.

Press ENTER to delete the data.

Press ENTER to confirm the deletion of data or press CLEAR to exit the function.

“REMOTE CONTROL”

The REMOTE CONTROL function is used to connect the SAS to an optional infrared data downloading device.

“CALIBRATION”

The “Calibration” function is a password protected menu accessible only by an authorised technician.

“AUTOSWITCH”

In order to save battery power, after 5 minutes of inactivity the instrument turns itself off.

With the “Autoswitch” function is possible to disable this feature.


“MULTI MODE” function

This programme is very useful for extending the time of sampling with the purpose of obtaining a more representative environmental sample “in actual operating conditions” (very important, for example, during surgery).


The total air volume to be sampled is aspirated with two or more sub volume aspirations (e.g. 1000 litres in ten runs of 100 litres at five minute intervals). Before entering the MULTI MODE programme you should, therefore, decide (a) total volume of air to be sampled onto the contact plate; (b) number of runs; (c) interval time.

Switch on the unit by touching the screen for at least two seconds.

Press MENU and select MULTI MODE function by pressing the “UP” or “DOWN” arrows.

Press ENTER to enter the multimode function.

Select the desired interval time between sampling by pressing the “UP” or “DOWN” arrows. Press the ENTER button to go to the next digit.

The interval time can be set from 0 to 59 minutes.

Press ENTER to confirm.

NUMBERS OF RUN			
[01] RUN			
^			
CLEAR	↓	↑	ENTER

Select the desired number of sampling cycles by pressing the “UP” or “DOWN” arrows.

Press the ENTER button to go to the next digit.

The maximum numbers of cycle is 19.

Press ENTER to confirm.

SINGLE RUN VOLUME			
[0000] Litres			
^			
CLEAR	↓	↑	ENTER

Select the desired volume of aspirated air for each sampling by pressing the “UP” or “DOWN” arrows.

Press the ENTER button to go to the next digit.

The maximum aspirated volume of the cycle is 1999 litres.

Press ENTER to confirm.

The display will show the TOTAL volume of air to be aspirated at the end of the cycles.


Check if it is the correct volume.

If not, press the CLEAR button and start again.


For example if you set 5 cycles of 200 litres each, you will aspirate a total volume of 1000 litres of air.

TOTAL RUN	
[1000] Litres	
CLEAR	ENTER


Pressing ENTER to confirm will set the SAS in multi mode.

Press START to run multi mode.

At the end of the run the SAS will display a flashing delay before the next cycle; number of completed cycles and residual cycles (runs) are also displayed.


“CONFIGURATION” function

To start with one of the eight operator programmable configurations memorised.


If the operator wants to use a pre-set configuration, (composed of site, identity, numbers of runs, numbers of cycles) up to 8 additional configurations can be programmed into the unit.

Switch on the unit by touching the screen for at least two seconds.


Press MENU and select CONFIGURATION function by pressing the “UP” or “DOWN” arrows.

Press ENTER to go to the function.

Select the SELECT function by pressing the "UP" or "DOWN" arrows and press ENTER.

Select the pre-set configuration desired (between 1 and 8) by pressing the "UP" or "DOWN" arrows and press ENTER to confirm.

Press START to run the sampling cycle.

Modify one of the eight operator programmable configuration memorised.

If you want to modify the pre-set configuration (composed of site, identity, number of runs, number of cycles) you must use the MODIFY function.

Switch the unit on by touching the screen for at least two seconds.

Press MENU and select CONFIGURATION function by pressing the "UP" or "DOWN" arrows.

Press ENTER to go to this function.

CONFIGURATION			
MODIFY			
CLEAR	↓	↑	ENTER

Select the SET function by pressing the “UP” or “DOWN” arrows and press ENTER.

MODIFY			
CONFIGURATION [1]			
CLEAR	↓	↑	ENTER

You must choose the configuration that you want to modify function by pressing the “UP” or “DOWN” arrows and press ENTER.

You can now modify the site name.

SITE			
[ABCDEFGH]			
	^		
CLEAR	↓	↑	ENTER

Change the letter indicated by “^” by pressing the “UP” or “DOWN” arrows and press ENTER to go to the next letter.

SITE			
[LABOEF GH]			
		^	
CLEAR	↓	↑	ENTER

Repeat until you reach the last letter and then confirm by pressing ENTER; you can now modify the identity.

IDENTIFY			
[12345678]			
^			
CLEAR	↓	↑	ENTER

Change the figure indicated by “^” by pressing the “UP” or “DOWN” arrows and press ENTER to go to the next figure.

IDENTIFY			
[35345678]			
^			
CLEAR	↓	↑	ENTER

Repeat until you reach the last figure and then confirm by pressing ENTER; you can now choose to programme the SAS in single mode or multi mode.

INTERVAL TIME			
[00] Minutes			
^			
CLEAR	↓	↑	ENTER

Select the desired interval time between sampling by pressing the “UP” or “DOWN” arrows. Press ENTER to go to the next digit.

The interval time can be set from 0 to 59 minutes.

Note: If you want to programme a single cycle set INTERVAL TIME to 0 minutes

Press ENTER to confirm.

NUMBER OF RUNS			
[01] RUN			
^			
CLEAR	↓	↑	ENTER

Select the desired numbers of sampling cycles by pressing the “UP” or “DOWN” arrows.

Press the ENTER button to go to the next digit.

The maximum numbers of cycle is 19.

Note: If you want to programme a single cycle set NUMBER OF RUNS to 1.

Press ENTER to confirm.

SINGLE RUN VOLUME			
[0000] Litres			
^			
CLEAR	↓	↑	ENTER

Select the desired volume of aspirated air for each sampling by pressing the “UP” or “DOWN” arrows.

Press the ENTER button to go to the next digit.

The maximum aspirated volume of cycle is 1999 litres.

Press ENTER to confirm.

The display will show the TOTAL volume of air to be aspirated at the end of the cycles.

Check if it is correct volume.

If not, touch the CLEAR button and start again.

For example if you set 5 cycles of 200 litres each, you will aspirate 1000 litres of air.

TOTAL RUN RECORD MODE	
[1000] Litres	
CLEAR	ENTER

Press ENTER to confirm and save the configuration.

“POWER OFF” function

Press MENU and select POWER OFF function by pressing the “UP” or “DOWN” arrows.


MENU			
POWER OFF			
CLEAR	↓	↑	ENTER

Press ENTER to switch the instrument off.

Note: The instrument will automatically switch off after 5 minutes.

Operations

Adjustable plate holders

The SAS SUPER ISO 100 and SAS SUPER ISO 180 air samplers are designed for use with standard contact or Petri dishes.

The plate holders are adjustable to allow for different brands of plates to be fitted. Adjustable plate holders in the sampler head allow the operator to use contact plates from various sources.

Troubleshooting

Review the information in the table below to troubleshoot operating problems.

Problem	Cause	Solution
Unit does not start	Instrument is not turned on	Touch the screen for at least two seconds
	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and if it needs replacing order Cat. No. 710-0973
	Battery is too old	Check battery and replace it
Battery discharges after few minutes working	Battery is low	Recharge battery
	Battery charger is not working	Check battery charger and if it needs replacing order Cat. No. 710-0973
	Battery is too old	Check battery and replace it
“LOW BATTERY” message	Battery is low	Recharge battery
“CALIBRATION EXPIRED” message	Instrument needs to be calibrated	Send the instrument to VWR International or an authorised dealer
Infrared remote switch doesn't switch the unit on	Remote control battery is low	Replace infrared remote battery
Microbiological media is dehydrated after sampling	Media is damaged	Control expiration date of media and check that the agar is not dehydrated before sampling
	Sampling time is too long	Shorten the sampling time

Accessories

Description	Cat. No.
Aspirating heads for contact plates, Ø 55 mm	
Stainless steel	710-0880
Aluminium	710-0892
Sterile daily	710-0890
Aspirating heads for Petri dishes, Ø 90 mm	
Stainless steel	710-0878
Aluminium	710-0886

Sterile daily head	710-0891
Other accessories	
Soft carrying case	710-0896
Aluminium carrying case	710-0875
Bio-Transport autoclavable carrying case	113-8185
Handle for Bio-Transport carrying case	113-8186
Floor tripod	710-0889
SAS-Holder table and wall stainless steel	710-0963
Battery charger with universal plug for both models	710-0973
Adapter* to convert contact plate model to accept 90 mm Petri dishes	710-0882
SAS stainless steel Petri head + adaptor	710-0877
SAS aluminium Petri head + adaptor	710-0879
IQ OQ PQ validation protocols for SAS Super ISO 100 and 180	710-0956
SAS software for downloading data from SAS Super ISO (to use with interface)	710-0970
Interface for SAS software for SAS Super ISO	710-0971
Infrared remote control for SAS Super ISO	710-0969

Technical services

Web resources

Visit the VWR website at www.vwr.com for:

- Complete technical service contact information
- Access to the VWR online catalogue and information about accessories and related products
- Additional product information and special offers

Contact us: For information or technical assistance contact your local VWR representative or visit.

www.vwr.com.

Warranty

VWR International warrants that this product will be free from defects in material and workmanship for a period of two (2) years from date of delivery. If a defect is present, VWR will, at its discretion and cost, repair, replace, or refund the purchase price of this product to the customer, provided it is returned during the warranty period. This warranty does not apply if the product has been damaged by accident, abuse, misuse, or misapplication, or from ordinary wear and tear. If the required maintenance and inspection services are not performed according to the manuals and any local regulations, such warranty will become invalid, except to the extent, the defect of the product is not due to such non performance.

Items being returned must be insured by the customer against possible damage or loss. This warranty shall be limited to the aforementioned remedies. IT IS EXPRESSLY AGREED THAT THIS WARRANTY WILL BE IN LIEU OF ALL WARRANTIES OF FITNESS AND IN LIEU OF THE WARRANTY OF MERCHANTABILITY.

Compliance with local laws and regulations

The customer is responsible for applying for and obtaining the necessary regulatory approvals or other authorisations necessary to run or use the product in its local environment. VWR will not be held liable for any related omission, or for not obtaining the required approval or authorisation, unless any refusal is due to a defect of the product.

Equipment disposal

This equipment is marked with the crossed out wheeled bin symbol to indicate that this equipment must not be disposed of with unsorted waste.

Instead it's your responsibility to correctly dispose of your equipment at lifecycle end by handing it over to an authorised facility for separate collection and recycling. It's also your responsibility to decontaminate the equipment in case of biological, chemical and/or radiological contamination, so as to protect from health hazards the persons involved in the disposal and recycling of the equipment.

For more information about where you can drop off your equipment, please contact your local dealer from whom you originally purchased this equipment.

By doing so, you will help to conserve natural and environmental resources and you will ensure that your equipment is recycled in a manner that protects human health.

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


The screen displays 'UTILITY MODE' at the top, followed by a horizontal line and the word 'SITE' in the center. Below 'SITE' are four buttons: 'CLEAR', a downward arrow, an upward arrow, and 'ENTER'.

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The screen displays 'SITE' at the top, followed by a horizontal line and the text '[ABCDEFGH]' in the center. Below the text is a small '^' symbol. Below the '^' are four buttons: 'CLEAR', a downward arrow, an upward arrow, and 'ENTER'.

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The screen displays 'SITE' at the top, followed by a horizontal line and the text '[LABOEF GH]' in the center. Below the text is a small '^' symbol. Below the '^' are four buttons: 'CLEAR', a downward arrow, an upward arrow, and 'ENTER'.

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A rectangular screen with a black border. At the top, the text "UTILITY MODE" is centered above a horizontal line. Below the line, the word "SITE" is centered. At the bottom, there are four rounded rectangular buttons: "CLEAR", a downward arrow, an upward arrow, and "ENTER".

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A rectangular screen with a black border. At the top, the word "SITE" is centered above a horizontal line. Below the line, the text "[ABCDEFGH]" is centered, with a small caret (^) symbol centered below it. At the bottom, there are four rounded rectangular buttons: "CLEAR", a downward arrow, an upward arrow, and "ENTER".

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A rectangular screen with a black border. At the top, the word "SITE" is centered above a horizontal line. Below the line, the text "[LABOEF GH]" is centered, with a small caret (^) symbol centered below it. At the bottom, there are four rounded rectangular buttons: "CLEAR", a downward arrow, an upward arrow, and "ENTER".

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A rectangular screen with a black border. At the top, the text "UTILITY MODE" is centered above a horizontal line. Below the line, the word "IDENTIFY" is centered. At the bottom, there are four rounded rectangular buttons: "CLEAR", a downward arrow, an upward arrow, and "ENTER".

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The screenshot shows a menu titled "SITE" with a horizontal line below it. Below the line, the text "[LABOEF GH]" is displayed, with a small caret (^) under the 'L'. At the bottom, there are four buttons: "CLEAR", a "DOWN" arrow, an "UP" arrow, and "ENTER".

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The screenshot shows a menu titled "IDENTIFY" with a horizontal line below it. Below the line, the text "[12345678]" is displayed, with a small caret (^) under the '1'. At the bottom, there are four buttons: "CLEAR", a "DOWN" arrow, an "UP" arrow, and "ENTER".

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The screenshot shows the same "IDENTIFY" menu, but now the text is "[35345678]" with the caret (^) under the '3'.

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The screenshot shows a menu titled "INTERVAL TIME" with a horizontal line below it. Below the line, the text "[00] Minutes" is displayed, with a small caret (^) under the first '0'. At the bottom, there are four buttons: "CLEAR", a "DOWN" arrow, an "UP" arrow, and "ENTER".

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VWR International AB
VWR International AG
Teknolojileri Ltd.Şti.
VWR International Ltd