

Datasheet

GPRIN2 mouse monoclonal antibody (hybridoma)

Catalog Number: H00009721-M

Regulation Status: For research use only (RUO)

Product Description: Mouse monoclonal antibody raised against a full-length recombinant GPRIN2.

Immunogen: GPRIN2 (AAH11672.1, 1 a.a. ~ 461 a.a) full-length recombinant protein with GST tag. MW of the GST tag alone is 26 KDa.

Sequence:

MSSSHPEPGPWAPLSPRLQPLSQSSSSLLGEGREQR
PELHKASSTMWQAQLGEASTRPQAPEEEGNPPESM
KPARASGPKARPSAGGHWRSSTVGNVSTMGGDLCL
RLRAPSAAMQRSHSDLVRSTQMRGHSGARKASLSC
SALGSSPVHRAQLQPGGTSGQGGQAPAGLERDLAPE
DETSNSAWMLGASQLSVPLDLGDTTAHSSSAQAEP
KAAEQLATTTCHALPPAALLCGMREMREVGAGGCCH
ALPATGILAFPKLVASVSESLQAQHGKIHCRSLSGGL
PGHSHCCAHLWGPAGLVPEPGSRTKDVWTMTSANDL
APAEASPLSAQDAGVQAAPVAACKAVATSPSLEAPAA
LHVFPEVTLGSSLEEAPSPVRDVRWDAEGMTWEVYG
AAVDPEVLGVAIQKHLEMQFEQLQRAPASEDSLSVEG
RRGPLRAVMQSLRRPSCCGCSGAAPE

Host: Mouse

Reactivity: Human

Applications: ELISA, WB-Re, WB-Tr
(See our web site product page for detailed applications information)

Protocols: See our web site at
<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

Entrez GeneID: 9721

Gene Symbol: GPRIN2

Gene Alias: GRIN2, KIAA0514, MGC15171