

## Datasheet

### GPR150 polyclonal antibody

**Catalog Number:** PAB25999

**Regulation Status:** For research use only (RUO)

**Product Description:** Rabbit polyclonal antibody raised against synthetic peptide of GPR150.

**Immunogen:** A synthetic peptide corresponding to 16 amino acids at C-terminus of human GPR150.

**Host:** Rabbit

**Reactivity:** Human

**Applications:** IHC-P

(See our web site product page for detailed applications information)

**Protocols:** See our web site at

<http://www.abnova.com/support/protocols.asp> or product page for detailed protocols

**Specificity:** BLAST analysis of the peptide immunogen showed no homology with other human proteins, except ATP2B2 (56%).

**Form:** Liquid

**Purification:** Immunoaffinity chromatography

**Recommend Usage:** Immunohistochemistry

(Formalin/PFA-fixed paraffin-embedded sections) (1 ug/mL)

The optimal working dilution should be determined by the end user.

**Storage Buffer:** In PBS (0.1% sodium azide)

**Storage Instruction:** Store at 4°C. For long term storage store at -80°C.

Aliquot to avoid repeated freezing and thawing.

**Entrez GeneID:** 285601

**Gene Symbol:** GPR150

**Gene Alias:** PGR11