

Name: Rabbit polyclonal anti-ZNF95 antibody
Product Data Sheet - ANTIBODY**Catalog: TA315268**

Components:	• Rabbit polyclonal anti-ZNF95 antibody (TA315268)
Amount:	100ul
Immunogen:	The antiserum was produced against synthesized peptide derived from internal of human ZNF95.
Host:	Rabbit
Isotype:	IgG
Species Reactivity:	Human
Guaranteed Applications:	IF
Suggested Dilutions:	IF: 1:100~1:500, ELISA: 1:5000
Concentration:	1mg/ml
Buffer:	Phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Storage Condition:	Shipped at -20C. Upon delivery store at -20C. Dilute in PBS (pH7.3) if necessary. Stable for 12 months from date of receipt. Avoid repeated freeze-thaws.


Target

Target Name:	Homo sapiens zinc finger with KRAB and SCAN domains 5 (ZKSCAN5), transcript variant 1
Alternative Name:	ZFP95; ZNF914; ZSCAN37
Database Link:	NP_055384
Function:	

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

© 2014 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850

Validation Data



Immunofluorescence analysis of A549 cells, using ZNF95 antibody. The picture on the right is treated with the synthesized peptide.

* More validation images may be available on our website: <http://www.origene.com/antibody/TA315268.aspx>

This product is to be used for laboratory only. Not for diagnostic or therapeutic use.

© 2014 OriGene Technologies, Inc. 9620 Medical Center Dr., Suite 200, Rockville, MD 20850
