



Material	26185.297
Material description	n-Pentane
Grade	AnalaR NORMAPUR
CAS Number	109-66-0
Molecular formula	C ₅ H ₁₂
Molecular mass	72.15

Characteristics	Specifications
Assay (on anhydrous substance)	Min. 99.0 %
Assay (C5 isomers)	Min. 99.7 %
Acidity	Max. 0.004 meq/g
Boiling point	35 --> 37 °C
Colouration	Max. 10 APHA
Density (20/4)	0.625 --> 0.632
Aromatics (as C ₆ H ₆)	Max. 50 ppm
Evaporation residue	Max. 10 ppm
Total S (as SO ₄)	Max. 50 ppm
Water	Max. 100 ppm

Signature

We certify that this batch conforms to the specifications listed above.

This document has been produced electronically and is valid without a signature.

Isabelle Habay, Head of Laboratory - Briare
VWR International S.A.S.; Z.I. de Vaugereau; FR-45250 Briare;
France

For Professional use in Laboratory or Manufacturing. Not for use as an Active Pharmaceutical Ingredient or Food or Animal Feed. Suitability and intended use of the product remains the responsibility of the user.