



VWR® Carboys

PP Carboys

Use for storing and dispensing solutions and media. Ideal for sterile water storage, polypropylene (PP) carboys can be autoclaved before refilling to stop recurring bacterial growth. Excellent chemical resistance.

HDPE Carboys

The high-density polyethylene (HDPE), disposable carboys are ideal for storing and transfer of sterile fluids and pharmaceutical/biotech reagents.

- Good for Cold storage (Up to -80° C), working temperature +120° C to -80° C
- High Impact strength even at negative temperature.
- Widely used for long storage due to High chemical resistance property
- Widely used in Food & Milk Industries
- Sterile Carboy available for high purity application.



LDPE Carboys

- Preferred for Hi purity application as it is high resistance to chemicals (Example: Widely used in Ophthalmic product range, like Eye droppers, nasal spray etc)
- More preferred resin for Cold storage (Up to -80° C), working temperature +80° C to -80° C
- Resistance to physical impacts (products will not break even if it falls down by accidentally, where as PP Carboys may lead to crack at negative temperature)
- LDPE Product will bounce back to its shape as its High flexible material.
- Widely used in Pharma industries.

Carboys have both wide mouth and narrow mouth option to suite the application
Graduated to identify the volume dispensed



- **Carboy:** Refer to table for material type, 100 % Virgin, FDA 21 CFR Compliance, USP Class VI passed.
- **Closure:** Polypropylene, 100% Virgin, FDA 21 CFR Compliance, USP Class VI passed with food-grade color.
- **BSE/TSE Statement:** These products are manufactured without the use of raw materials of animal origin.
- **Latex Free Statement:** No materials containing latex or natural rubber are used in the manufacturing, handling and packaging processes for these products.

VWR® Carboys

VWR® Rectangular Carboy with Spigot, Wide Mouth						
Material	Capacity	Cap Size	Neck ID	Sterility	Cat. No.	UOM
HDPE	10 L (2.6 gal.)	100 mm	85 mm	Nonsterile	216-1719	Each
PP	10 L (2.6 gal.)	100 mm	85 mm	Nonsterile	216-1735	Each
HDPE	20 L (5.3 gal.)	100 mm	85 mm	Nonsterile	216-1700	Each
PP	20 L (5.3 gal.)	100 mm	85 mm	Nonsterile	216-1732	Each
VWR® Round Carboy, Wide Mouth, Without Spigot						
LDPE	10 L (2.6 gal.)	100 mm	85 mm	Nonsterile	216-1708	Each
PP	10 L (2.6 gal.)	100 mm	85 mm	Nonsterile	216-1730	Each
PP	20 L (5.3 gal.)	100 mm	85 mm	Nonsterile	216-1733	Each
VWR® Round Carboy, Narrow Mouth, Without Spigot						
HDPE	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1736	Each
LDPE	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1704	Each
PP	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1703	Each
HDPE	20 L (5.3 gal.)	83 mm	65 mm	Sterile	216-1729	Each
HDPE	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1737	Each
LDPE	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1706	Each
PP	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1701	Each
LDPE	50 L (13.2 gal.)	83 mm	65 mm	Nonsterile	216-1738	Each
PP	50 L (13.2 gal.)	83 mm	65 mm	Nonsterile	216-1702	Each
VWR® Heavy-Duty Round Carboy, Narrow Mouth, Without Spigot						
PP	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1709	Each
PP	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1707	Each
VWR® Round Carboy, Narrow Mouth, With Spigot						
LDPE	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1731	Each
PP	10 L (2.6 gal.)	83 mm	65 mm	Nonsterile	216-1710	Each
LDPE	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1717	Each
PP	20 L (5.3 gal.)	83 mm	65 mm	Nonsterile	216-1705	Each
LDPE	50 L (13.2 gal.)	83 mm	65 mm	Nonsterile	216-1734	Each
PP	50 L (13.2 gal.)	83 mm	65 mm	Nonsterile	216-1718	Each
Accessories						
Stopcock for VWR Carboys					216-1745	Each

Caution

1. Before Autoclaving separate the caps and Spigot / Stopcock threads.
2. Do not stack products, this may cause vacuum and may lead to physical damage/collapse of product.